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Sommaire

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
FHN Daily Trials

HD Conventiionelle 3x4h versus HDDQ 6x2h

Level I : lowering UFR ± lowering blood flow
 Level II : administration of saline solution but not lowering UFR
 Level III : administration of saline solution + lowering UFR

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Values are expressed as n (%), mean ± SD. Abbreviations: LVM, left ventricular mass, LVEF, left ventricular ejection fraction.

LVM: Left Ventricular Ejection Fraction
 LVEF: Left Ventricular Ejection Fraction

1. FHN Trial Group. In-Center Hemodialysis. N Engl J Med. 2010;363:2287-300.
 3. Fessi H et al. Physidia S³: RECAP Study. J Clin Med. 2023;12:1357.
 4. FHN Trial Group. Frequent HD & BP. Hemodial Int. 2015;19(3):386-401.

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

Figure 2 | Minutes to recovery after a hemodialysis session in the Daily Trial. F/up, follow-up.


2. Garg AX et al. Frequent Hemodialysis QoL. Kidney Int. 2017;91:746-754.
 3. Fessi H et al. Physidia S3: RECAP Study. J Clin Med. 2023;12:1357.

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
- ▶ Citation dans Nutrients 2017:

"On intensive dialysis schedules, higher efficiency and/or frequency may make an unrestricted diet possible (dialysis liberalizes diet)"⁶

- ▶ Citation dans Hemodialysis International 2024:

"More than 65% of HHD patients stated that they were not bothered at all by limitations on food and water intake"⁵

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


Figure 2 — Survival among home hemodialysis patients, matched in-center hemodialysis patients, and unmatched in-center hemodialysis patients, in intention-to-treat analysis following peritoneal dialysis technique failure.


7. McCullough PA et al. Intensive HD & LVH. Am J Kidney Dis. 2016;68(Suppl 1)
 8. Weinhandl ED et al. Home vs In-Center HD Survival. J Am Soc Nephrol. 2012;23:895-904.
 9. Kansal SK et al. Home vs In-Center HD Transplant. Perit Dial Int. 2019;39:25-34.

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Differences in hospitalization rates across dialysis modalities and probability of a Monday/ Tuesday admission

Kidney360

		Dialysis modality	Overall rate of admissions/ 1000 patient days	Odds of experiencing a Monday/Tuesday admission
<p>Canadian Organ Replacement Register</p> <p>Maintenance Dialysis (Canada)</p> <p>Hospitalizations (were attributed to the previous treatment if they occurred within 30 days of a treatment change)</p>	<p>111748 Hospitalization episodes</p> <p>n=27430 individuals</p> <p>Differences in the proportion of patients experiencing a Monday/ Tuesday admission with all other days of the week were compared</p>	<p>In-center Hemodialysis (ICHD) n=21,095</p>	2.98	<i>Reference</i>
		<p>Frequent In-center Hemodialysis n=392</p>	3.76	OR 0.97 (95% CI 0.89-1.07)
		<p>Home Hemodialysis (HHD) n=236</p>	2.16	OR 0.99 (95% CI 0.88-1.12)
		<p>Frequent Home Hemodialysis n= 445</p>	2.13	OR 0.89 (95% CI 0.81-0.97)
		<p>Peritoneal Dialysis n=5,262</p>	2.71	OR 0.95 (95% CI 0.93-0.97)
		<p>The OR was lower when restricted to hospitalization episodes for CV reasons comparing frequent HHD to ICHD (OR 0.68, 95% CI 0.48-0.96)</p>		

Conclusions In this nationally representative cohort, we identified that the probability of a Monday/Tuesday admission was lower for frequent HHD and PD compared with ICHD, most notably for frequent HHD and hospitalizations due to cardiovascular causes. Gaining a better understanding of the reasons behind this observation may help to develop future strategies to reduce overall and cause specific hospitalization for dialysis patients.

Karthik K. Tennankore, Annie-Claire Nadeau-Fredette, Kara Matheson, et al. **Home Versus In-center Dialysis and Day of the Week Hospitalization: A Cohort Study.** *Kidney360*. 10.34067/KID.0003552021
Visual Abstract by Edgar Lerma, MD, FASN

7. McCullough PA et al. Intensive HD & LVH. *Am J Kidney Dis.* 2016;68(Suppl 1)

10. Tennankore K et al. Home vs In-Center HD & Hospitalization. *Kidney360*. 2022;3(1):103-12

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
11. Walker RC, Howard K, Morton RL. Home hemodialysis review: patient-centered & economic. Clinicoecon Outcomes Res. 2017;9:149-161.

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
“ There were no significant differences in the RR for urea, phosphate, kFLC and albumin, whereas the weekly RRs for creatinine and $\beta 2$ -microglobulin were significantly higher with the Physidia S³ system than those obtained with OL-HDF¹⁵ ”

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Références

1. FHN Trial Group. In-Center Hemodialysis: Six Times per Week versus Three Times per Week. *N Engl J Med.* 2010;363:2287-300.
2. Garg AX, et al. Patients receiving frequent hemodialysis have better health-related quality of life compared to patients receiving conventional hemodialysis. *Kidney Int.* 2017;91:746-754.
3. Fessi H, Nicoud P, Serrato T, Gilbert O, Courivaud C, Daoud S, et al. Two Years' Experience of Intensive Home Hemodialysis with the Physidia S³ System: Results from the RECAP Study. *J Clin Med.* 2023;12:1357.
4. FHN Trial Group. Effects of Frequent Hemodialysis on Blood Pressure: Results from the Randomized Frequent Hemodialysis Network Trial. *Hemodial Int.* 2015;19(3):386-401.
5. Watanabe Y, Tsubakihara Y, Nishi S, Akiba T, Suzuki M, Ohira S, et al. Home hemodialysis and conventional in-center hemodialysis in Japan: A comparison of health-related quality of life. *Hemodial Int.* 2014;18(Suppl 3)
6. Piccoli GB, Sofroniadou S, Trippoli S, Capizzi I, Cabiddu G, Maxia S, et al. The diet and haemodialysis dyad: Three eras, four open questions and four paradoxes. *Nutrients.* 2017;9(4):372.
7. McCullough PA, Chan CT, Depner TA, Levin NW, Greene T, Beck GJ, et al. Intensive hemodialysis, left ventricular hypertrophy, and cardiovascular disease. *Am J Kidney Dis.* 2016;68(Suppl 1)
8. Weinhandl ED, Gilbertson DT, Collins AJ. Survival in daily home hemodialysis and matched thrice-weekly in-center hemodialysis patients. *J Am Soc Nephrol.* 2012;23:895-904.
9. Kansal SK, Trentini A, MacRae JM, Davenport A, Blake PG. Survival and kidney transplant incidence on home versus in-center hemodialysis, following peritoneal dialysis technique failure. *Perit Dial Int.* 2019;39:25-34.
10. Tennankore K, Soroka SD, Chan CT. Home versus in-center dialysis and day of the week hospitalization: A cohort study. *Kidney360.* 2022;3(1):103-12.
11. Walker RC, Howard K, Morton RL. Home hemodialysis review: patient-centered & economic. *Clinicoecon Outcomes Res.* 2017;9:149-161.
12. Fessi H, Szelag JC, Courivaud C, Daoud S, Dillinger JG, Canaud B, et al. Safety and efficacy of short daily hemodialysis with Physidia S³ system: Clinical performance assessment during the training period. *J Clin Med.* 2022;11(8):2123.
13. Potier J, Queffeuilou G, Bouet J. Are all dialyzers compatible with the convective volumes suggested for postdilution online hemodiafiltration? *Int J Artif Organs.* 2016;39(9):460-70.
14. Thomas M, Vincent E, Todorova V. P1051A Single Session and Weekly Beta 2 Microglobulin Removal During Short Frequent Hemodialysis Using Physidia Technology. *Nephrol Dial Transplant.* 2020;35(Suppl_3).
15. Moya MG, Kammer M, Canaud B, Fleuret J, de Seigneux S, Thomas M. Elimination of middle and large uremic toxins in short daily home hemodialysis with low dialysate volume: A randomized crossover clinical trial. *Nephrol Dial Transplant.* 2022;37(Suppl_3).
16. Frequent Hemodialysis Network Trial Group. In-center hemodialysis and blood pressure outcomes. *J Am Soc Nephrol.* 2006;17(2):546-55.

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